

Gibson Assembly® Site-Directed Mutagenesis Kit

Quick Reference Manual for Experienced Users

For additional product information, an online SDM primer design tool, and detailed protocol visit sgidna.com/sdm-kit

Catalog Numbers GA2100-S and GA2100-10

Products and Storage Conditions

Gibson Assembly® Site-Directed Mutagenesis Kit Contents

Quantity	Component	Cat. GA2100-S (5 Reactions)	Cat. GA2100-10 (10 Reactions)	Storage Temperature
		Volume		
1 tube (each)	GA SDM PCR Mix (2X)	315 µL	625 µL	-20°C
	GA SDM Assembly Mix A (2X)*	25 µL	50 µL	
	GA SDM Assembly Mix B (2X)*	50 µL	100 µL	
	Control Plasmid (5 ng/µL)	5 µL (2 Control Rxns)	5 µL (2 Control Rxns)	
	Control Primer Mix A (25 µM)	5 µL	5 µL	
	Control Primer Mix B (25 µM)	5 µL	5 µL	

*Aliquot the Assembly Mixes to reduce the number of freeze-thaw cycles.

Mutagenesis Amplification Reaction Setup

1. Thaw GA SDM PCR Mix (2X) on ice.
2. For each fragment, prepare the following reaction mixture in a thin-walled PCR tube on ice:

Component	Volume
GA SDM PCR Mix (2X)	12.5 µL
dsDNA template (100 pg – 30 ng)	X µL
SDM Forward Primer (5 µM)	2.5 µL
SDM Reverse Primer (5 µM)	2.5 µL
PCR water	7.5–X µL
Total	to 25 µL

3. (Optional) For the positive control, prepare the following reaction mixtures in thin-walled PCR tubes. Set up one reaction for Control Amplicon A and set up a separate reaction for Control Amplicon B (i.e. use Control Primer Mix A in one reaction and Control Primer Mix B in a separate reaction.)

Component	Volume
GA SDM PCR Mix (2X)	12.5 µL
Control Plasmid (5 ng/ µL)	1.0 µL
Control Primer Mix A or B (25 µM)	1.0 µL
PCR water	10.5 µL
Total	25 µL

4. Mix by pipetting up and down.

Mutagenesis Amplification Reaction

1. Transfer tubes to a thermocycler. Program and start the thermocycler using the following conditions:

Denature	98°C	1 min	1 Cycle
Amplify	98°C	10–30 sec	30 Cycles
	Primer T _m *	30 sec	
	72°C	30 sec per kb	
Extend	72°C	5 min	1 Cycle
Hold	4°C	as necessary	1 Cycle

*Control Primers T_m is 60°C.

2. After the reaction is complete, analyze 5 µL of the product on a 0.8% agarose gel.
3. Minimize original template carryover by:
 - Gel-purifying the PCR product
 - or
 - Treating the PCR product with DpnI and column-purifying the DpnI-treated reaction.

Gibson Assembly® Reaction Setup

1. Thaw GA SDM Assembly Mix A (2X) on ice.
2. Adjust the concentration of the purified, mutagenized fragments according to the following table. In a PCR tube, combine the insert(s) and vector, bringing the total volume to 5 µL with nuclease-free water.

Fragment	Size	Amount	Insert:Vector Molar Ratio
Insert	≤ 1 kb	10–40 ng	1:1 to 5:1
	1–5.5 kb	10–25 ng	1:1 to 2:1
Vector	2–10 kb	25 ng	—

Note: For dilute fragments, you may increase the total combined sample volume to 10 µL. For 10 µL samples, increase subsequent assembly volumes 2-fold.

3. (Optional) For the positive control, combine 25 ng of gel- or column-purified Positive Control Amplicon A with 25 ng of gel- or column-purified Positive Control Amplicon B. Bring the total volume to 5 µL with nuclease-free water.
4. Vortex GA SDM Assembly Mix A (2X) immediately before use (after it is thawed).
5. In a tube on ice, combine 5 µL of volume-adjusted, purified mutagenized fragments and 5 µL of thawed GA SDM Assembly Mix A (2X).
6. Mix the reactions by pipetting up and down.
7. Vortex and spin down all reactions. Proceed to the Gibson Assembly® Reaction.

The Gibson Assembly® Reaction

1. Transfer tubes to a thermocycler. Program and start the thermocycler using the following conditions (1 cycle):

3' end chew back	37°C for 5 minutes
Inactivation	75°C for 20 minutes
Slow cooling	0.1°C/second to 60°C
Annealing	60°C for 30 minutes
Slow cooling	0.1°C/second to 4°C

2. Thaw GA SDM Assembly Mix B (2X) on ice. **Vortex and briefly centrifuge the thawed master mix immediately before use.**
 3. While keeping the tubes on ice, add 10 µL of GA SDM Assembly Mix B (2X) to the completed reaction from step 1.
 4. Mix the reaction by pipetting.
 5. Incubate samples using the following conditions:
- | | | |
|--------|---------|---------------------|
| Repair | 1 cycle | 45°C for 15 minutes |
|--------|---------|---------------------|
6. Store completed reactions at –20°C or proceed to downstream applications (e.g. transformation).
 7. (Optional) Electrophorese 50% of the assembly reaction (e.g. 10 µL of a 20 µL reaction) on a 0.8–2% agarose gel. A high molecular weight smear is indicative of a successful assembly reaction.

Transformation Guidelines

For electrocompetent cells

- Use 1–2.5 µL of Gibson Assembly® Reaction per 30–50 µL of competent cells.
- Use high efficiency competent cells (transformation efficiency $\geq 1 \times 10^9$ CFU/µg pUC19). We recommend TransforMax™ EPI300™ Electrocompetent *E. coli* (Epicentre® Cat. No. EC300110)

For chemically competent cells

- Use 1 µL of Gibson Assembly® Reaction per 50 µL of competent cells.
- Use high efficiency competent cells (transformation efficiency $\geq 1 \times 10^9$ CFU/µg pUC19). We recommend *E. coli* 10G Chemically Competent Cells (Lucigen Cat. No. 60107)

Plating Guidelines

Following transformation, plate 1/2–1/50 of the transformation reaction onto an LB plate containing appropriate antibiotics. We recommend plating 2 volumes. For the positive control, plate 1/100 volume of the transformed reaction onto an LB plate with 100 µg/mL ampicillin or carbenicillin and 40 µg/mL X-Gal and 0.1 mM IPTG.

Additional product information and resources are available at sgidna.com/sdm-kit

Technical Services techservices@sgidna.com

Trademarks *E. coli*® is a registered trademark of Lucigen Corporation.
TransforMax™ and EPI300™ are trademarks of Epicentre Technologies Corporation.
Gibson Assembly® is a registered trademark of SGI-DNA, Inc.

Label License
Research Use Only Except as otherwise agreed in writing by our authorized representative, this product is for INTERNAL RESEARCH USE ONLY AND NOT FOR HUMAN, ANIMAL, THERAPEUTIC OR DIAGNOSTIC USE. For additional information about your rights under this research license, please see our website at sgidna.com.

Disclaimer The material in this manual is for informational purposes only and is subject to change without prior notice at any time. SGI-DNA and/or its affiliates assume no responsibility for any errors that may appear in this document.

Gibson Assembly® US Patent Nos. 7,776,532, 8,435,736, and 8,968,999.